

# Product Data Sheet

<b>Product name</b>	Rollator
<b>REF</b>	RL200
<b>General Information</b>	The rollator features four wheels, a padded seat, and a storage bag located under the seat. It also comes with ergonomic hand grips with brakes for added safety and stability. Its lightweight and foldable design make it convenient for mobility and storage.
<b>Specifications</b>	Light Weight Aluminum - 7.5kg
	Foam padded seat
	Practical storage bag
	Lightweight and foldable design
	Four wheels for stability
	Ergonomic hand grips with manual brakes
	Dimensions: 54,7 x 62,7 cm
	User weight: 136 kg

## Product photo



### Headquarters:

Medlogics.com LTD  
Lyuben Karavelov 6  
9000 Varna, Bulgaria

**E** info@medlogics.com  
**W** www.medlogics.com

 Medlogics.com